



MKRSIII[®]

Sistema Modular Para Revision de Rodilla



Movimiento como símbolo de vida



Componente Femoral

ANVISA 10417940046



CoCrMo

Direito
Ref. Nº

04.15.08.00001
04.15.08.00002
04.15.08.00003
04.15.08.00004

Esquerdo
Ref. Nº


04.15.09.00001
04.15.09.00002
04.15.09.00003
04.15.09.00004

Tam.

P
M
G
EG

A x B

64 x 54mm
69 x 58mm
72 x 62mm
77 x 66mm



Platô Tibial

ANVISA 10417940041



UHMPWE

Ref. N°	Base	Fêmur	Ref. N°	Base	Fêmur	Ref. N°	Base	Fêmur	Tam.
04.17.06.01010	P	P / M	04.17.07.02010	M	G	04.17.06.03010	G	G / EG	10mm
04.17.06.01012	P	P / M	04.17.07.02012	M	G	04.17.06.03012	G	G / EG	12mm
04.17.06.01015	P	P / M	04.17.07.02015	M	G	04.17.06.03015	G	G / EG	15mm
04.17.06.01018	P	P / M	04.17.07.02018	M	G	04.17.06.03018	G	G / EG	18mm
04.17.06.01021	P	P / M	04.17.07.02021	M	G	04.17.06.03021	G	G / EG	21mm
04.17.06.01025	P	P / M	04.17.07.02025	M	G	04.17.06.03025	G	G / EG	25mm
04.17.06.02010	M	P / M	04.17.07.03010	G	M	04.17.06.04010	EG	G / EG	10mm
04.17.06.02012	M	P / M	04.17.07.03012	G	M	04.17.06.04012	EG	G / EG	12mm
04.17.06.02015	M	P / M	04.17.07.03015	G	M	04.17.06.04015	EG	G / EG	15mm
04.17.06.02018	M	P / M	04.17.07.03018	G	M	04.17.06.04018	EG	G / EG	18mm
04.17.06.02021	M	P / M	04.17.07.03021	G	M	04.17.06.04021	EG	G / EG	21mm
04.17.06.02025	M	P / M	04.17.07.03025	G	M	04.17.06.04025	EG	G / EG	25mm


Base Tibial

ANVISA 10417940061



CoCrMo

Ref. N°	Tam.	A x B
04.14.07.00001	P	67 x 42mm
04.14.07.00002	M	67 x 46mm
04.14.07.00003	G	74 x 46mm
04.14.07.00004	EG	74 x 50mm



Pino Extensor

ANVISA 10417940061



CoCrMo

Reto Ref. N°	Ø x Compr.
04.25.02.11070	11 x 70mm
04.25.02.12070	12 x 70mm
04.25.02.13070	13 x 70mm
04.25.02.14070	14 x 70mm
04.25.02.15070	15 x 70mm
04.25.02.16070	16 x 70mm
04.25.02.11100	11 x 100mm
04.25.02.12100	12 x 100mm
04.25.02.13100	13 x 100mm
04.25.02.14100	14 x 100mm
04.25.02.15100	15 x 100mm
04.25.02.16100	16 x 100mm



CoCrMo

Angulado Ref. N°	Ø x Compr.
04.25.03.11070	11 x 70mm
04.25.03.12070	12 x 70mm
04.25.03.13070	13 x 70mm
04.25.03.14070	14 x 70mm
04.25.03.15070	15 x 70mm
04.25.03.16070	16 x 70mm
04.25.03.11100	11 x 100mm
04.25.03.12100	12 x 100mm
04.25.03.13100	13 x 100mm
04.25.03.14100	14 x 100mm
04.25.03.15100	15 x 100mm
04.25.03.16100	16 x 100mm

Calço Femoral

ANVISA 10417940046



CoCrMo

Distal Ref. N°	Tam.	Altura
04.05.12.05001	P	05mm
04.05.12.05002	M	05mm
04.05.12.05003	G	05mm
04.05.12.05004	EG	05mm
04.05.12.10001	P	10mm
04.05.12.10002	M	10mm
04.05.12.10003	G	10mm
04.05.12.10004	EG	10mm



CoCrMo

Posterior Ref. N°	Tam.	Altura
04.05.11.0500	P	05mm
1	M	05mm
04.05.11.0500	G	05mm
2	EG	05mm
04.05.11.0500	P	10mm
3	M	10mm
04.05.11.0500	G	10mm
4	EG	10mm
04.05.11.1000		
1		
04.05.11.1000		
2		
04.05.11.1000		



CoCrMo

Anterior Ref. N°	Tam.
04.05.10.00001	P
04.05.10.00002	M
04.05.10.00003	G
04.05.10.00004	EG


Trava para fixação
no componente
femoral

Calço Tibial

ANVISA 10417940061



CoCrMo

1/3 Angulado Ref. N°	Tam.	Ângulo
04.05.09.22001	P	22
04.05.09.22002	M	°
04.05.09.22003	G	22
04.05.09.22004	EG	°
		22
		°
		22
		°



CoCrMo

Total Angulado Ref. N°	Tam.	Ângulo
04.05.08.07001	P	7
04.05.08.07002	M	°
04.05.08.07003	G	7
04.05.08.07004	EG	°
		7
		°
		7
		°



CoCrMo

Parcial Angulado Ref. N°	Tam.	Ângulo
04.05.06.16001	P	16°
04.05.06.16002	M	16°
04.05.06.16003	G	16°
04.05.06.16004	EG	16°
04.05.06.26001	P	26°
04.05.06.26002	M	26°
04.05.06.26003	G	26°
04.05.06.26004	EG	26°



CoCrMo

Parcial Paralelo Ref. N°	Tam.	Altura
04.05.07.05001	P	05mm
04.05.07.05002	M	05mm
04.05.07.05003	G	05mm
04.05.07.05004	EG	05mm
04.05.07.10001	P	10mm
04.05.07.10002	M	10mm
04.05.07.10003	G	10mm
04.05.07.10004	EG	10mm


Trava para fixação
na base tibial